

Construction & Performance of an Indexed Lateral Flow Microarray (iLFM™) Strip



Proteins were printed onto CN-95 nitrocellulose using the Lighthouse Diagnostics Staccato piezo dispensing system. Spots were developed with gold anti-Hu IgG4 conjugate under lateral flow conditions. Consistent spot morphology was achieved across 1, 20, and 50 drop deposition of 427 pL/drop on 450 μ m pitch. The strip comprises 5x6 sub-arrays for dietary food and environmental allergens.



Mean Signal Intensity: 50 drops, 89.1 ± 5.6 (6.3%); 20 drops, 70.6 ± 4.4 (6.3%); 1 drop, 19.8 ± 2.3 (11.6%).



For more information: contact_us@quantiscientifics.com

Staccato Bass SDC 90-B




Piezo dispensing technology



427 pL droplets



High positioning accuracy



5 x 6 indexed array




427 pL droplets



450 μ m pitch

NUMBERED SPOT MAP – DIAMETER (μ m)



20 drops


50 drops

LEGEND



SUMMARY

Total Spots	Mean Diameter (μ m)	Range (μ m)	Std Dev (μ m)	CV (%)	Edge-to-Edge Spacing (μ m)
24	340	325 – 350	9.7	2.8	~110



Sample
(Au⁴⁰ anti-IgG4)

Mean Spot Dia. 340 μ m, CV ~ 2.8 %

All diameters in micrometers (μ m)
Estimated $\pm 20 \mu$ m

SMALL SPOTS – DIAMETER MAP (μ m)

All diameters in micrometers (μ m) • Estimates $\pm 15 \mu$ m



1 drop

LEGEND



SMALL SPOT SUMMARY

Total Spots	Mean Diameter (μ m)	Range (μ m)	Std Dev (μ m)	CV (%)	Edge-to-Edge Spacing* (μ m)
24	89	82 – 99	4.6	5.2	~360



* Edge-to-edge spacing = Pitch (450 μ m) – Mean Diameter

Mean Spot Dia. 89 μ m CV ~ 5.2 %

The small spots mirror the large spot pattern on the same 450 μ m pitch grid.